

Product Datasheet

CHEK1 antibody (orb107026)

Description	Rabbit polyclonal to CHEK1. Serine/threonine-protein kinase which is required for checkpoint-mediate...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CHEK1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	54434
Uniprot ID	O14757
NCBI	NP_001107593.1 , NP_001231775.1 , NP_001107594.1 , NP_001265.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

