



Product Datasheet

BLVRA antibody (orb107029)

Description	Rabbit polyclonal to BLVRA. Reduces the gamma-methene bridge of the open tetrapyrrole, biliverdin IX...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	BLVRA
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	33428
Uniprot ID	P53004
NCBI	NP_000703.2, NP_001240752.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 cell line
...