



Product Datasheet

PITPNA antibody (orb107030)

Description	Rabbit polyclonal to PITPNA. Catalyzes the transfer of PtdIns and phosphatidylcholine between membra...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	PITPNA
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	31806
Uniprot ID	Q00169
NCBI	NP_006215.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...