



Product Datasheet

ZNF434 antibody (orb107040)

Description	Rabbit polyclonal to ZNF434. May be involved in transcriptional regulation.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic Peptide
Target	ZSCAN32
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	78728
Uniprot ID	Q9NX65
NCBI	NP_060280.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of MCF-7 cel...