

Product Datasheet

LILRA3 antibody (orb107041)

Description	Rabbit polyclonal to LILRA3. May act as soluble receptor for class I MHC antigens.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	LILRA3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	47472
Uniprot ID	Q8N6C8
NCBI	NP_006856.3 , NP_001166125.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

