


Product Datasheet

FAM159A antibody (orb107042)

Description	Rabbit polyclonal to FAM159A. The exact function of this protein remains unknown.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic Peptide
Target	SHISAL2A
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	20306
Uniprot ID	Q6UWV7
NCBI	NP_001036158.1
Dilution Range	IF/ICC: 1:25, WB: 1:1000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...

