

Product Datasheet

CABLES2 antibody (orb107045)

Description	Rabbit polyclonal to CABLES2. Unknown. Probably involved in G1-S cell cycle transition.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CABLES2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	52235
Uniprot ID	Q9BTV7
NCBI	NP_112492.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

