

Product Datasheet

TSPYL2 antibody (orb107053)

Description	Rabbit polyclonal to TSPYL2. Part of the CASK/TRB1/TSPYL2 transcriptional complex which modulates ge...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TSPYL2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	79435
Uniprot ID	Q9H2G4
NCBI	NP_071400.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

