

Product Datasheet

ZNF71 antibody (orb107069)

Description	Rabbit polyclonal to ZNF71. ZNF71 may be involved in transcriptional regulation.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ZNF71
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	54498
Uniprot ID	Q9NQZ8
NCBI	NP_067039.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

