

Product Datasheet

TRIM24 antibody (orb107071)

Description	Rabbit polyclonal to TRIM24. The protein encoded by this gene mediates transcriptional control by in...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TRIM24
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	116831
Uniprot ID	O15164
NCBI	NP_056989.2 , NP_003843.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

