

Product Datasheet

STAT1 antibody (orb107084)

Description	Rabbit polyclonal to STAT1. The protein encoded by this gene is a member of the STAT protein family....
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	STAT1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	87335
Uniprot ID	P42224
NCBI	NP_009330.1 , NP_644671.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

