



Product Datasheet

DMRTA2 antibody (orb107085)

Description	Rabbit polyclonal to DMRTA2. DMRTA2 may be involved in sexual development.
Species/Host	Rabbit
Reactivity	Human, Mouse, Zebrafish
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic Peptide
Target	DMRTA2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	53356
Uniprot ID	Q965C8
NCBI	NP_115486.1
Dilution Range	WB: 1:1000, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.



Confocal immunofluorescence analysis of ...