



Product Datasheet

CNBP antibody (orb107086)

Description	Rabbit polyclonal to CNBP. This gene encodes a nucleic-acid binding protein with seven zinc-finger d...
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	CNBP
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	19463
Uniprot ID	P62633
NCBI	NP_001120668.1, NP_001120665.1, NP_003409.1, NP_001120666.1, NP_001120664.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...