



Product Datasheet

YBX1 antibody (orb107089)

Description	Rabbit polyclonal to YBX1. Mediates pre-mRNA alternative splicing regulation. Binds to splice sites ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic Peptide
Target	YBX1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	35924
Uniprot ID	P67809
NCBI	NP_004550.2
Dilution Range	WB: 1:1000, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.



Confocal immunofluorescence analysis of ...