


Product Datasheet

GFP antibody (orb10709)

Description	GFP antibody
Species/Host	Rabbit
Reactivity	All
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	KLH conjugated synthetic peptide derived from belt jellyfish Green Fluorescent protein (25-120/238 aa)
Target	GFP
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P42212
Dilution Range	WB=1:1000-10000, IF=1:100-500, IHC-P=1:100-500, IHC-F=1:100-500
Expiration Date	12 months from date of receipt.



Paraffin-embedded pig liver tissue fixed...