

Product Datasheet

POU4F2 antibody (orb107090)

Description	Rabbit polyclonal to POU4F2. POU4F2 is a member of the POU-domain family of transcription factors. P...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	POU4F2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	43087
Uniprot ID	Q12837
NCBI	NP_004566.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

