


Product Datasheet

BASP1 antibody (orb1071)

Description	BASP1 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr
Immunogen	KLH conjugated synthetic peptide derived from human BASP1/Nap22 (3-100/227 aa)
Target	BASP1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	25 kDa
Uniprot ID	P80723
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...