



Product Datasheet

MLXIPL antibody pair (orb107110)


Description	Rabbit polyclonal to MLXIPL. MLX-interacting protein-like (MLXIPL) is a nuclear protein involved in ...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ChIP, IP, WB
Immunogen	-
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	95
Uniprot ID	Q99MZ3
Expiration Date	12 months from date of receipt.



Western blot (WB) of total cell extracts...



c) human HeLa, (b)



d) human Jurkat, (c) human HepG2; using ...