



Product Datasheet

TRMT2A antibody (orb107112)

Description	Rabbit polyclonal to TRMT2A. TRMT2A function is unknown. However, TRMT2A is orthologous to the mouse...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	N-terminal region of human TRMT2A
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	69
Uniprot ID	Q8IZ69
Expiration Date	12 months from date of receipt.



Western blot of total cell extracts from...



IP-WB.
Immunoprecipitation (IP) was perf...