



Product Datasheet

TICAM1 antibody (orb107116)

Description	Rabbit polyclonal to TICAM1. TIR domain-containing adapter molecule 1 (TICAM1) is Involved in innate...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human TICAM1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	96
Uniprot ID	Q8IUC6
Expiration Date	12 months from date of receipt.



Western blot (WB)
of total cell
extracts...



IP-WB.
Immunoprecipitation
(IP) was perf...