



Product Datasheet

PPP2R2C antibody (orb107119)

Description	Rabbit polyclonal to PPP2R2C. Serine/threonine-protein phosphatase 2A 55 kDa regulatory subunit B ga...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	C-terminal region of human PPP2R2C
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	52
Uniprot ID	Q9Y2T4
Expiration Date	12 months from date of receipt.



Western blot of total cell extracts from...



IP-WB.
Immunoprecipitation (IP) was perf...