


Product Datasheet

Rat Endothelin 1/EDN1 ELISA Kit (orb107175)

Description	Rat Endothelin 1/EDN1 ELISA Kit
Reactivity	Rat
Range	3.9 pg/ml - 250 pg/ml
Tested Applications	ELISA
Immunogen	A peptide reflecting the target protein in human. Sequence: C53-W73
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	23135 MW
Uniprot ID	P22388
Sensitivity	0.5 pg/ml
Expiration Date	Please enquire.



Rat
Endothelin
PicoKine
ELISA Kit
standa...