



Product Datasheet

CD116 antibody (orb10732)


Description	CD116 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from mouse GM-CSFR alpha (101-230/388 aa)
Target	GM-CSFR alpha
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	38 kDa
Uniprot ID	P15509
Dilution Range	IHC-P=1:100-500, IF=1:100-500, IHC-F=1:100-500
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of human cancerous bone marrow cells ...



IF of human cancerous bone marrow cells ...



Paraffin-embedded mouse intestine tissue...