



Product Datasheet

DUSP4 Antibody (orb107384)

Description	DUSP4 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q13115
Expiration Date	12 months from date of receipt.



Western blot analysis
of Gel: 10%SDS-
PAGE...



Immunohistochemical
analysis of formalin...