



Product Datasheet

RAB3c Antibody (orb107462)

Description	RAB3c Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q96E17
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse lung tiss...



Immunohistochemical analysis of formalin...