



Product Datasheet

CRMP5 Antibody (orb107488)

Description	CRMP5 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9BPU6
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of Gel: 10%SDS-PAG...