



Product Datasheet

LAIR2 Antibody (orb107541)

Description	LAIR2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q6ISS4
Expiration Date	12 months from date of receipt.



Western blot analysis
of Gel: 10%SDS-
PAG...



Immunohistochemical
analysis of formalin...