



Product Datasheet

LPAR4 Antibody (orb107567)

Description	LPAR4 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q99677
Expiration Date	12 months from date of receipt.



Western blot analysis
of Human
leiomyosa...



Immunohistochemical
analysis of formalin...