



Product Datasheet

Map2 Antibody (orb107587)

Description	Map2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P11137
Expiration Date	12 months from date of receipt.



Western blot analysis
of Gel:
8+10+12%SDS-PAGE



Immunohistochemical
analysis of formalin-fixed paraffin-embedded (FFPE) tissue