


Product Datasheet

MCF2L antibody (orb107644)

Description	Goat polyclonal to MCF2L
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	MCF2L
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 150kDa band observed in lysates of cell line Jurkat and of Human Thymus (calculated MW of 131kDa according to NP_001106203.1). Recommended concentration: 1-3µg/ml.
Protein Sequence	QLDNQATVQR
Clonality	Polyclonal
MW	131; 127
Entrez	23263
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat (A) and Human Thymus (B) lysates.