


Product Datasheet

MMP14 antibody (orb107657)


Description	Goat polyclonal antibody to MMP14
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	MMP14
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 60kDa band observed in lysates of cell lines NIH3T3, A549 and HeLa (calculated MW of 65.9kDa according to NP_004986.1). Recommended concentration: 1-3µg/ml. Primary incubation was 1 hour.. IHC: Paraffin embedded Human Liver. Recommended concentration: 5µg/ml.
Protein Sequence	ESATPLRFREVP
Clonality	Polyclonal
MW	65.9
Entrez	4323
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 lysate (35...



Western blot analysis of NIH3T3, (Lane 1...



Immunohistochemical staining of Human Li...