


Product Datasheet

CCDC3 antibody (orb107666)

Description	Goat polyclonal to CCDC3
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	CCDC3 (aa134-148)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: In transfected HEK293 transiently expressing Human CCDC3 (DYKDDDDK tagged), a band of approx. 36kDa is observed. No bands are observed in mock-transfected HEK293 and the same band is observed using anti-DYKDDDDK antibody. Recommended concentration: 0.5-1. IHC: Paraffin embedded Human Colon. Recommended concentration: 5µg/ml.
Protein Sequence	PHGVNFQDAIFPDTQ
Clonality	Polyclonal
MW	30.7
Entrez	83643
Dilution Range	ELISA: 1:64000, WB: 0.5-1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293 overexpr...