

Product Datasheet

Netrin-G1 antibody (orb107671)

Description	Goat polyclonal to Netrin-G1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA
Target	Netrin-G1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Preliminary experiments gave an approx 29kDa band in Human, Mouse and Rat Brain lysates from different regions after 0.1µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the c
Protein Sequence	SLYGQLDTTKLRD
Clonality	Polyclonal
MW	60.5; 53.9; 49.3
Entrez	22854
Dilution Range	ELISA: 1:128000, WB: 0.1 µg/ml
Expiration Date	12 months from date of receipt.