


Product Datasheet

LAT2 antibody (orb107674)

Description	Goat polyclonal to LAT2
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	LAT2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 26kDa band observed in Human Peripheral Blood Lymphocytes lysates (calculated MW of 26.6kDa according to NP_054865.2). Recommended concentration: 0.01-0.03µg/ml.
Protein Sequence	RSEKIYQQRSLRED
Clonality	Polyclonal
MW	26.6
Entrez	7462
Dilution Range	ELISA: 1:128000, WB: 0.01-0.03 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Peripheral Blood Lymphocytes lysates.