

Product Datasheet

TRIP4 antibody (orb107676)

Description	Goat polyclonal to TRIP4
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Target	TRIP4 / ASC1 (aa94-106)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Preliminary experiments in Human Testis, Thyroid gland and Ovary lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with o
Protein Sequence	DGQKSGDHLKRGR
Clonality	Polyclonal
MW	66.1
Entrez	9325
Dilution Range	ELISA: 1:128000, WB: 1 µg/ml
Expiration Date	12 months from date of receipt.