


Product Datasheet

IFNAR1 antibody (orb107680)

Description	Goat polyclonal to IFNAR1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	IFNAR1 (aa353-364)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 65kDa band observed in lysates of cell line K562 (calculated MW of 63.5kDa according to NP_000620.2). Recommended concentration: 0.3-1µg/ml.
Protein Sequence	PKQSGNTPVIQD
Clonality	Polyclonal
MW	63.5
Entrez	3454
Dilution Range	ELISA: 1:128000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 lysate (35...)