

Product Datasheet

IFNGR1 antibody (orb107681)

Description	Goat polyclonal to IFNGR1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	IFNGR1 (aa181-193)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 70kDa band observed in lysates of cell line HepG2 (calculated MW of 54.4kDa according to NP_000407.1). The observed molecular weight corresponds to the glycosylated form. Recommended concentration: 0.3-1µg/ml.
Protein Sequence	SEIQYKILTQKED
Clonality	Polyclonal
MW	54.4
Entrez	3459
Dilution Range	ELISA: 1:128000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.

