


Product Datasheet

Human Amphiregulin protein (orb107683)

Description	Amphiregulin is an EGF related growth factor and was originally isolated from the conditioned media ...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	Amphiregulin
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	SVRVEQVVKP PQNKTESENT SDKPKRKKKG GKNGKNRRNR KKKNPCNAEFQNFCHGECK YIEHLEAVTC KCQQEYFGER CGEKSMKTHS MIDSSLK
Purity	> 95% by SDS-PAGE and HPLC analyses
MW	11.3kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
Amphiregulin
...