


Product Datasheet

Human 4-1BB Ligand protein (orb107685)

Description	Recombinant Human 4-1BBL is a member of the tumor necrosis factor (TNF) receptor family. This recept...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	4-1BB Ligand
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	REGPELSPDD PAGLLDLRQG MFAQLVAQNV LLIDGPLSWY SDPGLAGVSLTGGLSYKEDT KELVAKAGV YYVFFQLELR RVVAGEGSGS VSLALHLQPLRSAAGAAALA LTVDLPPASS EARNSAFGFQ GRLLHLSAGQ RLGVHLHTEARARHAWQLTQ GATVLGLFRV TPEIPAGLPS PRSE
Purity	> 95 % by SDS-PAGE and HPLC analyses.
MW	19.4 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human 4-
1BB Ligand

...