

Product Datasheet

Murine BD-1 protein (orb107687)

Description	Defensin are 3-4 kDa antimicrobial peptides of which three distinct families have been identified; a...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	BD-1
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	VGILTSLGRR TDQYKCLQHG GFCLRSSCPN NTKLQGTCKP DKPNCKKS
Purity	> 97% by SDS-PAGE and HPLC analyses
MW	5.2 kDa
Expiration Date	6 months from date of receipt.