

Product Datasheet

Rat BD-1 protein (orb107688)

Description	Defensin are 3-4 kDa antimicrobial peptides of which three distinct families have been identified; D...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	BD-1
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in 20mM PB, 500mM NaCl, pH 7.0.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	DQYRCLQNGG FCLRSSCP SH TKLQGTCKPD KPNCCRS
Purity	> 95% by SDS-PAGE and HPLC analyses
MW	4.1 kDa
Expiration Date	6 months from date of receipt.