


Product Datasheet

Rat BD-4 protein (orb107690)

Description	Defensins (alpha and beta) are cationic peptides with a broad spectrum of antimicrobial activity tha...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	BD-4
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in 10mM PB, pH 7.4, 500mM NaCl.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	QSINNPITCL TKGGVCWGPC TGGFRQIGTC GLPRVRCCKK K
Purity	> 95% by SDS-PAGE and HPLC analyses
MW	4.4KDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Rat BD-4 protein