

Product Datasheet

Rat C-10 protein (orb107693)

Description	Chemokine (C-C motif) ligand 6 (CCL6) is a small cytokine belonging to the CC chemokine family that ...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	C-10
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in 20mM PB, pH 7.4, 50mM NaCl.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	GLIQDTVKED RPFNPTIIHQ GFQDSSDCCF SYASQIPCSR FIYYFPTSGG CTKPGIIFVTRKRKRVCANP SDQRVQTCIS TLKLGPRSGN SAIA
Purity	> 97% by SDS-PAGE and HPLC analyses
MW	10.4 kDa
Expiration Date	6 months from date of receipt.