

Product Datasheet

Rat CINC-2 alpha protein (orb107696)

Description	The three GRO cDNAs encode 107 amino acid precursor proteins from which the N-terminal 34 amino acid...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	CINC-2 alpha
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in 20mM PB, pH 7.4, 50mM NaCl.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	RELRCQCLKT LPRVDFENIQ SLTVTPPGPH CTQTEVIATL KDGQEVCLNP QAPRLQKIIQKLLKSDKSS
Purity	> 95% by SDS-PAGE and HPLC analyses.
MW	7.8 kDa
Expiration Date	6 months from date of receipt.