


Product Datasheet

Human CYR61 protein (orb107703)

Description	Cysteine-rich angiogenic inducer 61 (Cyr61) encoded by the Cyr61 gene is a dynamically expressed, mu...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	CYR61
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in citrate buffer solution, 300 mM NaCl, pH3.0.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	TCPAACHCPL EAPKCAPGVG LVRDGC GCCK VCAKQLNEDC SKTQPCDHTKGLECNFGASS TALKGICRAQ SEGRPCEYNS RIYQNGESFQ PNCKHQCTCI DGAVGCIPLCPQELSLPNLG CPNPRLVKVT GQCCEEWVCD EDSIKDPMED QDGLLGKELG FDASEVELTRNNELIAVGKG SSLKRLPVFG MEPRILYNPL QGQKCIVQTT SWSQCSKTCG TGISTRVTNDNPECRLVKET RICEVRPCGQ PVYSSLKKGK KCSKTKKSPE PVRFTYAGCL SVKKYRPHYKCGSCVDGRCCT PQLTRTVKMR FRCEDGETFS KNVMMIQSCK CNYNCPHANEAAFPFYRLFN DIHKFRD
Purity	> 95 % by SDS-PAGE and HPLC analyses
MW	39.4 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
CYR61
protein