

Product Datasheet

Human ENA-78 protein (orb107705)

Description	Epithelial cell-derived neutrophil-activating peptide 78 (ENA-78) is a member of the CXC subfamily o...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	ENA-78
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	LRELRCVCLQ TTQGVHPKMI SNLQVFAIGP QCSKVEVVAS LKNGKEICLD PEAPFLKKV IQKILDGGNKEN
Purity	> 95% by SDS-PAGE and HPLC analyses
MW	7.8 kDa
Expiration Date	6 months from date of receipt.