


Product Datasheet

Recombinant Murine Eotaxin protein (orb107706)

Description	Eotaxin also called CCL11 is a CC chemokine that signals through the CCR3 receptor. It is produced b...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	Eotaxin
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	HPGSIPTSCC FIMTSKKIPN TLLKSYKRIT NNRCTLKAIV FKTRLGKEIC ADPKKKWVQDATKHLDQKLQ TPKP
Purity	> 95% by SDS-PAGE and HPLC analyses
MW	8.4 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis
of
Recombinant Murine
E...