


Product Datasheet

Rat Eotaxin protein (orb107707)

Description	Eotaxin/CCL11 is a potent eosinophil chemoattractant that was originally purified from bronchoalveol...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	Eotaxin
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in 20mM PB, pH 7.4, 150mM NaCl.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	HPGSIPTSCC FTMTSKKIPN TLLKSYKRIT NNRCTLKAIV FKTKLGKEIC ADPKKKWVQDATKHLDQKLQ TPKP
Purity	> 95% by SDS-PAGE and HPLC analyses
MW	8.4 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Rat Eotaxin
protein