


Product Datasheet

Murine FGF-acidic protein (orb107711)

Description	Recombinant Murine Acidic Fibroblast Growth Factor. Murine aFGF, encoded by the FGF1 gene, is a memb...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	FGF-acidic
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	FNLPLGNYKK PKLLYCSNGG HFLRILPDGT VDGTRDRSDQ HIQLQLSAESAGEVYIKGTE TGQYLAMDTE GLLYGSQTPN EECLFLERLE ENHYNTYTSKKHAEKNWFGV LKKNNGSCKRG PRTHYGQKAI LFLPLPVSSD
Purity	> 96% by SDS-PAGE and HPLC analyses
MW	15.8 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Murine FGF-acidic p...