


Product Datasheet

FGF-4 protein (orb107712)

Description	FGF-4 is a heparin binding growth factor that is a member of the FGF family. Proteins of this family...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	FGF-4
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
FGF-4 protein